



Prometheus NT.Plex: Differential Scanning fluorimetry using Intrinsic Protein Fluorescence (NanoDSF): Capillary Chips

PR-AC002 Prometheus NT.Plex nanoDSF Grade Standard Capillary Chips

- Standard capillary chips for proteins concentrations of at least 200 $\mu\text{g}/\text{mL}$ (with tryptophan).
- 2 x 8 Prometheus nanoDSF Grade Standard 24-Capillary Chips.

PR-AC006 Prometheus NT.Plex nanoDSF Grade High Sensitivity Capillary Chips

- High sensitivity capillary chips for protein concentrations below 200 $\mu\text{g}/\text{mL}$.
- 2 x 8 Prometheus nanoDSF Grade High Sensitivity 24-Capillary Chips.



last edited: 2024-04-04